

ENT TREATMENT UNIT ACCESSORIES

Otopront Automatic Nasosinus Irrigation System

Automatic irrigation of the nose and sinuses from the treatment unit



The Automatic Nasosinus Irrigation System rinses the nasal cavity and paranasal sinuses directly from the Otopront ENT treatment unit. Its ergonomic handpiece with two nasal olives is placed at both nostrils, and the irrigation fluid is delivered automatically and evenly, so rinsing is comfortable for the patient and leaves the doctor's other hand free. It is used, for example, in the treatment of chronic rhinosinusitis and in follow-up care after nasal and sinus surgery.

Key features

- Ergonomic handpiece with two anatomically shaped nasal olives for both nostrils
- Automatic, even delivery of the irrigation fluid (e.g. saline solution) from the treatment unit
- Comfortable for the patient, efficient for the doctor: one-handed operation
- Stainless steel and black handle with holder clip; braided supply hose
- Designed for integration into Otopront ENT treatment units
- Configuration and technical details on request

Request a quote

WhatsApp / Phone: +995 597 33 77 62

E-mail: david.bakhtadze@futuretechnology.info

Future Technology Ltd. · 68 Giorgi Chubinashvili St., Chugureti District, Tbilisi 0102, Georgia

otopront.ge · futuretechnology.info



Scan for details