

ENT TREATMENT UNIT ACCESSORIES

Otopront Adapter for All Pharmacy Bottles

Spray medication straight from the original pharmacy bottle



With the adapter for pharmacy bottles, medication can be sprayed directly from its original pharmacy bottle using the spray system of the Otopront ENT treatment unit, with no decanting into separate spray bottles. Spray heads fit standard pharmacy bottles of different sizes and work with Otopront single-use spray cannulas, for hygienic, precisely dosed application.

Key features

- Fits standard pharmacy bottles of different sizes
- Spray medication straight from the original bottle, with no decanting
- Stainless steel spray heads and bottle holders
- Works with Otopront single-use spray cannulas (green for aqueous, white for oily/viscous solutions)
- For the compressed-air spray handpiece of the Otopront ENT treatment unit
- Configuration and technical details on request

Request a quote

WhatsApp / Phone: +995 597 33 77 62

E-mail: david.bakhtadze@futuretechnology.info

Future Technology Ltd. · 68 Giorgi Chubinashvili St., Chugureti District, Tbilisi 0102, Georgia

otopront.ge · futuretechnology.info



Scan for details