

ENT TREATMENT UNIT ACCESSORIES

Otopront Tonsillar Lavage Device

Conservative cleaning of infected tonsils



The Otopront tonsillar lavage device (Tonsil Cleaner) cleans infected tonsils as a conservative treatment option and can give patients fast relief from pain and infection. Pressurised water or saline is sprayed onto the tonsils to loosen the infected coating, which is then removed through a tonsil-shaped suction tip with adjustable vacuum. Clean tonsils also make it easier to apply medication.

Key features

- Conservative treatment option for chronically infected tonsils
- Spray of water or saline (NaCl) loosens the infected coating
- Tonsil-shaped suction tip removes the material with adjustable vacuum from the unit
- Connects to the water/NaCl reservoir and suction of the Otopront ENT treatment unit
- Stainless steel handpiece; the complete handpiece is autoclavable
- Makes subsequent application of medication easier

Request a quote

WhatsApp / Phone: +995 597 33 77 62

E-mail: david.bakhtadze@futuretechnology.info

Future Technology Ltd. · 68 Giorgi Chubinashvili St., Chugureti District, Tbilisi 0102, Georgia

otopront.ge · futuretechnology.info



Scan for details